

QESTLab Release Notes

(c) 2008-2011 Spectra Quality Engineering Software Tools P/L

Spectra QEST
Suite 7, 132 O'Connell Street
North Adelaide, SA 5006 Australia

Phone: +61 8 8360 8900
Facsimile: +61 8 8360 8988
Support: +61 8 8360 8999



Overview

From time to time Spectra QEST releases service packs for the QESTLab. These are required for the inclusion of new features and the resolution of software defects. It is always only necessary to install the most recent service pack as it will include all previous changes.

While every effort is made to ensure service packs are free from error Spectra QEST can not guarantee this to be the case. If possible it is recommended that the service pack be implemented in a test environment before being installed on the production system

Some features described in the release notes may be available subject to licensing arrangements and therefore not accessible to all users.

Release Notes for QESTLab Version 3.3.1400

CR#978 Invoice copy text not showing on Invoice

Coffey Information:Cof11.161

The Coffey custom invoice header has been updated to allow for additional room on the top to display the various labels (as configured in the QEST Administrator Console).

Kindly contact the Spectra QEST HelpDesk for the updated custom header.

Bug#1512 Maximum Dry Density [RTA T130]

Coffey Information:Cof11.146

The Additive Proportion can now be entered and is reported to the nearest 0.1 percent. In addition, the Oversize Sieve and Percentage are now reported.

Bug#1520 Previous Page icon not shown for Lists

Coffey Information:Cof11.162

The "Page Up" icon, for lists, now appears as required.